



Digital Cognitive
Rehabilitation Therapy

NEUROvitalis

Prevention and Therapy
in Geriatrics

Training Progress and Outcome of Therapy

How HeadApp / NEUROvitalis documents the progress of the therapy and proves the outcome through the training:




Progress

Then program **documents every step in the training** completely automatically. For each client it is saved, when and how often they exercised and what the results of the individual tasks were.

This data is also collected when the patient works in home training and is available everywhere.

In this way it is achieved that the patient continues on a new training day where he finished his task on the last day. The patient does not have to solve a task twice and is neither overburdened nor under-challenged.




Outcome

The **outcome** of the training is measured with the **integrated screening**. A status is determined before the training, after a training period the status is checked again.


Screening contains different tasks and challenges than training, so screening results do not show any learning effects, but a real improvement in performance.

The second option to measure the outcome is to use the integrated **IADL questionnaire**.


Training Progress




All training data is automatically collected and saved for each patient. The data from training at the facility and at home are merged. Results are displayed by clicking the **"Results"** button.



After clicking **"Results"** a **history** is shown in a calendar view. Each color box is a training session. Every color stand for a defined training task. Hovering a box with the mouse cursor shows an info about the training. Click the box to show more details.



Clicking then button **"Correct and level"** shows a bar graph. The color bars hold the information about the training module used and the performance of the client (correct solutions in %). The gray bars show the difficulty of the task used.



Clicking then button **"Duration"** shows a second bar graph. These bars holding the information of training duration (minutes) for every session.

All data can be printed and exported as CSV file (for use in a statistical evaluation).

Outcome of Therapy

Option 1: Screening with pre- and post- tests

When the therapist uses the screening as pre- test, he gets valuable information about the deficits of the patient. The screening suggests, which therapy modules to use.

On a second point in time (after some weeks of therapy) the screening can be used again (in a parallel version). Now the outcome of therapy is shown.

Example printout:

NEUROvitalis Screening results for Steve Thomson, born 1955

| Screening module | Result | Date | Assessment | Percentile rank | |
|---|---|------------|-------------------|-----------------|-------------------|
| Word Pairs 1 Verbal Short-term Memory | memory capacity short | 9/22/2021 | impaired | 3% | ← Test result one |
| | | 11/10/2021 | lower average | 8% | ← Test result two |
| Test date one → | eye-hand speed | 9/22/2021 | lower average | 13% | |
| | | 11/10/2021 | normal | 25% | |
| Test date two → | Pathfinder Test Visual-spatial Attn. | 9/22/2021 | normal | 17% | |
| | | 11/10/2021 | normal | 23% | |
| | laterality | 9/22/2021 | lower average | 13% | |
| | | 11/10/2021 | normal | 16% | |
| Image Builder Spatial Cognition | solution speed | 9/22/2021 | strongly impaired | 1% | |
| | | 11/10/2021 | impaired | 4% | |
| Number Sort Working Memory | memorized numbers | 9/22/2021 | lower average | 8% | |
| | | 11/10/2021 | lower average | 10% | |
| Wordfinding Speech | word finding | 9/22/2021 | lower average | 11% | |
| | | 11/10/2021 | normal | 18% | |
| Tower of London Planning | problem solving | 9/22/2021 | lower average | 8% | |
| | | 11/10/2021 | lower average | 12% | |
| Word Pairs 2 Verb. Medium-term Memory | memory capacity medium | 9/22/2021 | impaired | 6% | |
| | | 11/10/2021 | lower average | 10% | |

The picture shows, that the client reaches a higher percentile rank (good outcome) after six weeks of therapy.

Option 2: IADL questionnaire

The **Lawton Instrumental Activities of Daily Living (IADL) Scale** is part of HeadApp's client administration.

The screenshot displays the 'Patient Administration' interface in HeadApp. At the top left is the HeadApp logo. The title 'Patient Administration' is centered at the top. To the right of the title is a search bar with the placeholder text 'search from list' and two checkboxes: 'Return to login' and 'Data protection'. Below the title, the 'currently selected patient: John Smith' is shown. Underneath, there are 'sort by:' options: 'last therapy session' (selected), 'Name / Code', and 'Year of birth'. A list of patients is shown below, with 'John Smith; born 1955; last therapy at 11/11/2021' highlighted. To the right of the patient list are several action buttons: 'Continue', 'New', 'Add email', 'Edit', 'Delete', 'IADL' (highlighted with a red box and showing 'Score = 3'), and 'Manage license'. At the bottom left, there are 'Back' and 'Max.' buttons. At the bottom center, there is a 'Therapist / institution' button. Below the patient list, there is a summary: 'Patient(s): 1 total, 1 treated during the last 4 weeks, 0 at home. Your license is valid for 20 patients in therapy until the 12/2/2021.'

This IADL scale is an appropriate instrument to assess independent living skills (Lawton & Brody, 1969).

The instrument is most useful for identifying how a person is functioning at the present time, and to identify **improvement or deterioration over time**. There are eight domains of function measured with the Lawton IADL scale.

Clients are scored according to their highest level of functioning in every category. A summary score ranges from 0 (low function, dependent) to 8 (high function, independent).